

Welcome to the Virtual Platform site for

## iTelemed

[Go Back Home](#)

Already have an account? [Login here](#)



Video Visiting



SMS and E-mail Notifications



Telephone consultation



Screen sharing and presentation



Secure Document sharing



ICD Integration

Welcome to the Virtual Platform site for

# TELEMEDICINE iTelemed

>> Book Now <<

Already have an account? [Login here](#)

Patient Clicks to Book Appointment

Administrator Login  
Doctors Login  
Staff Login

**Main Screen** gives access for patients to Book appointments, for Administrators to Manage the site, for staff to manage specific Clinic or branch, for Doctors to manage their encounters and schedules.

# Branches/ Clinics Screen



Login

- Displays branches/ clinics in both Map and List selection Formats.

Patient can select the clinic/ Branch they need to book Their appointment at. Either By clicking on select or on the Respective pin on the map.

Map Satellite  
intelemed platform-1.pdf

London

**Huron**  
1472 Huron St.  
519-851-6445

**BMC Baseline**

Select


Select

Click on "Pin" or "Select" to choose Clinic/Branch for appointment

# Service Selection Screen



Login



Click on the required service  
e.g. General Practitioner, ENT,  
Specialized Doctors, etc...

How would you like to have your appointment?



- Select Type of Appointment
- In Person: Patient to visit the clinic.
- Phone: Patient attends Encounter with Dr. over The phone.
- Video: Patient attends Encounter remotely using MEDXPRESS Software of APP.



# Reason for Visit Screen




Login

## Appointment

Edit


Baseline  
/w Dr. Practitioner  
Video



Every step is Registered  
in this section where patient  
Can edit and make  
modifications

### Reason for visit

Medication Refill/Renew	COVID-19 Related	UTI
STI/STD	Ear Infection	Eye Infection/Pinkeye
Sore Throat	Strep Throat	Back Pain
Abdominal Pain	Chest Pain	Fever/Cough/Flu
Bloodwork Requisition	Rash/Derma/Skin Irritation	Diarrhea/Constipation
Nausea/Vomiting	Headache/Migraine	Dizziness/Vertigo
Reason not listed?	Other	



Back

Patient can select the reason for the visit  
Mostly listed in the selection screen above, if not  
Patient can simply type what are the symptoms  
Or other reasons for the visit.

Next



# Date of Visit Screen

## Appointment Edit

Baseline  
/w Dr. Practitioner  
Video (Chest Pain)

Every step is Registered  
in this section where patient  
Can edit and make  
modifications

### Date of visit

July 2022

Today

Mon	Tue	Wed	Thu	Fri	Sat	Sun
-	-	-	-	1	2	3
4	5	6	7	8	9	10
11	12	13	14	15	16	17
18	19	20	21	22	23	24
25	26 Available times: 0	27 Available times: 11	28 Available times: 11	29 Available times: 11	30 Available times: 0	31 Available times: 0

Back Next

Patient can select the date to book the Appointment, using above interactive calendar With color coding defining availability and possibility to book or not.

# Time of Visit Screen

## Appointment

Edit

Baseline  
/w Dr. Practitioner  
Video (Chest Pain)

Every step is Registered  
in this section where patient  
Can edit and make  
modifications

Select a date and time

2022  
Wed, July 27

Pre. Day Today Next Day

5PM - 5PM (0 openings)

5PM - 5PM (0 openings)

5PM - 5PM (0 openings)

5PM - 5PM (0 openings)

**5PM - 6PM (5 openings)**

6PM - 7PM (6 openings)

5:10 PM 5:20 PM 5:30 PM

5:40 PM 5:50 PM

Home Login

>

Sat	Sun
2	3
9	10
16	17
23	24
30 Available times: 0	31 Available times: 0

Patient can select time slot and time to book  
Appointment depending on availability.

Next



## Appointment

[Edit](#)

Baseline  
/w Dr. Practitioner  
Video (Chest Pain)  
Jul. 27/2022 - 5:40 PM

Every step is Registered  
in this section where patient  
Can edit and make  
modifications

### Already have an account?

Email or Phone Number

admin

Password

.....

[Forgot your password?](#)

Login

[Continue as new user](#)

Back

Patient are asked to log in, using their account if  
already registers, or to create a new account.


# Patient Info (account creation) Screen



Logout

## Appointment Edit

Baseline  
/w Dr. Practitioner  
Video (Chest Pain)  
Jul. 27/2022 - 5:40 PM



Every step is Registered  
in this section where patient  
Can edit and make  
modifications

### Patient Info

Choose Patient \*

Title: Please select


First Name \* [ ]

Last Name \* [ ]

Gender \* Please select

Date of Birth: 1990-01-01

← Back Next →



Patient use this screen to complete their info.

# Patient medical coverage Screen



Logout

## Appointment

Edit

Baseline  
/w Dr. Practitioner  
Video (Chest Pain)  
Jul. 27/2022 - 5:40 PM

## Contact Info

Edit

said-reda@esa-solutions.com  
519-864-8563

## Patient Info

Edit

Said Reda  
Dec. 31/1969  
Male

Every step is Registered  
in this section where patient  
Can edit and make  
modifications

## Health Card

I don't have a valid Health Card

I have a valid Health Card

Please Select Your Status

Canadian with no Health Card

International Student

Visitor or Out of Country

← Back

Next →

Patient use this screen to complete their medical coverage info, including national or provincial or private health insurance coverage or no coverage.



**Appointment** [Edit](#)

Baseline  
/w Dr. Practitioner  
Video (Chest Pain)  
Jul. 27/2022 - 5:40 PM

**Contact Info** [Edit](#)

said-reda@esa-solutions.com  
519-864-8563

**Patient Info** [Edit](#)

Said Reda  
Dec. 31/1969  
Male

### Health Card



I don't have a valid Health Card  I have a valid Health Card

Ontario

Health Card Number: 1213-134-442

Version Code: ee

[← Back](#) [Next →](#)



Every step is Registered in this section where patient can edit and make modifications

Patient use this screen to complete their medical coverage info, inputting data of their national or provincial or private health insurance coverage



**Appointment** Edit

Baseline  
/w Dr. Practitioner  
Video (Chest Pain)  
Jul. 27/2022 - 5:40 PM

**Contact Info** Edit

said-reda@esa-solutions.com  
519-864-8563

**Patient Info** Edit

Said Reda  
Dec. 31/1969  
Male

**Health Card** Edit

Ontario  
1213-134-442-ee

### Patient Address

Address1 \*  Unit # (If Applicable)

City \*  Postal Code \*  Province \*

← Back Next →

Patient use this screen to complete home address data.

Every step is Registered in this section where patient can edit and make modifications



**Appointment** [Edit](#)

Baseline  
/w Dr. Practitioner  
Video (Chest Pain)  
Jul. 27/2022 - 5:40 PM

**Contact Info** [Edit](#)

said-reda@esa-solutions.com  
519-864-8563

**Patient Info** [Edit](#)

Said Reda  
Dec. 31/1969  
Male

**Health Card** [Edit](#)

Ontario  
1213-134-442-ee

**Address** [Edit](#)

910 Roulston St

## Medical Info

Current Family Doctor (if applicable)

Dr. Smith

Any Medical Allergies?

Peniciline

Patient use this screen to provide info about family Doctor, allergies, if any.

← Back

Next →

Every step is Registered in this section where patient Can edit and make modifications



**Appointment** [Edit](#)

Baseline  
/w Dr. Practitioner  
Video (Chest Pain)  
Jul. 27/2022 - 5:40 PM

**Contact Info** [Edit](#)

said-reda@esa-solutions.com  
519-864-8563

**Patient Info** [Edit](#)

Said Reda  
Dec. 31/1969  
Male

**Health Card** [Edit](#)

Ontario  
1213-134-442-ee

**Address** [Edit](#)

910 Roulston St

## Pharmacy Info

Do you have a preferred Pharmacy?  Yes  No

Pharmacy of Preference \*

Shoppers

Pharmacy Address \*

1365 Huron Street London ON N5Y 6M7

Patient use this screen to provide info about pharmacy they prefer to use.

[← Back](#)

[Next →](#)

Every step is Registered in this section where patient can edit and make modifications



**Appointment** [Edit](#)

Baseline  
/w Dr. Practitioner  
Video (Chest Pain)  
Jul. 27/2022 - 5:40 PM

**Contact Info** [Edit](#)

said-reda@esa-solutions.com  
519-864-8563

**Patient Info** [Edit](#)

Said Reda  
Dec. 31/1969  
Male

**Health Card** [Edit](#)

Ontario  
1213-134-442-ee

**Address** [Edit](#)

910 Roulston St

## Confirm Appointment?

Comments

fever at night 38.5 C

I agree with [iTeled Med Terms & Conditions](#) \*

Patient use this screen, as a last step to confirm Appointment. Can add comment related to the Appointment to make doctor of staff aware of.

[← Back](#)

[Confirm](#)

Every step is Registered in this section where patient Can edit and make modifications



# Patient Scheduled appointment Screen

Upload File

Patient can upload files in any format.

New Booking

Patient can create new appointment

Upcoming

Past

Tuesday, Jul 26

2022

10:15am ● (Validating Health Card) Dr. Dr. for Said Reda



Wednesday, Jul 27

2022

5:40pm ● (Confirmed) Dr. Practitioner for Said Reda



This screen shows summary of scheduled appointments, past, present and future including time, status (color coded )

Patient can view their appointment type. (in Person, Phone Call, or Video Call)

Patient can reschedule appointment.

Patient can cancel appointment

# Patient files upload screen

**Add New File**

**General Information**

User or Client  
Said Reda (1969)

File(s) **Choose Files** No file chosen

« Back

Open

This PC > Pictures >

Search Pictures

Organize New folder

Said

This PC

Desktop

Documents

Downloads

MediaCenter-

Music

Pictures

Camera Roll

Camera Roll

Saved Pictures

Screenshots

pexels-tima-miroshnichenko-8376302.jpg

File name:

All files (\*.\*)

Open Cancel

Patient can upload files in any format.

# Patient Profile screen

## Clients



Show 50 entries

Search:

First Name	Last Name	Date of Birth	Health Card Number	Version Code	Please Select Your Status	Patient Phone Number	action
Said	Reda	31-12-1969	1213-134-442	ee	HC	519-864-8563	<a href="#">Edit</a>

Showing 1 to 1 of 1 entries

Previous 1 Next

Patient can view, edit and make changes to their profile/

Upload File

New Booking

Upcoming

Past

Wednesday, Jul 27

2022

5:40pm ● (Confirmed) Dr. Practitioner for Said Reda



Patient can start video encounter thru appointments screen, clicking on the video icon.

# System Administrator Home Screen

Search...

All Branches

All Doctors

Select Doctors

July 25 - 31, 2022

New Booking

Month

Week

3-Day

Day

Current Week

	Mon, Jul 25	Tue, Jul 26 (B:2, C:0)	Wed, Jul 27 (B:1, C:0)	Thu, Jul 28	Fri, Jul 29	Sat, Jul 30	Sun, Jul 31
8am							
8:05							
8:10							
8:15							
8:20							
8:25							
8:30							
8:35							
8:40							
8:45							
8:50							
8:55							
9am							
9:05							
9:10							
9:15							
9:20							
9:25							
9:30							
9:35							
9:40							
9:45							
9:50							
9:55							
10am							
10:05							
10:10							
10:15							
10:20							

Clinic schedule (color coded based on working days and hours)

Scheduled appointment (type: Video, Phone, In Person)  
Color coded based on Appointment status

# Schedule/ Booking List Screen

+ New Booking

Day

<< Back July 30, 2022 Next >>

Current Day

Show 50 entries Search:

↑↓	Status	Time	Name	Unpaid	Reason	Action	Assign	Edit
	Canceled	30-07-2022 10:00 AM to 10:15 AM	Said Reda (1969)	\$0	Back Pain	Select		Edit
	Confirmed	30-07-2022 10:20 AM to 10:35 AM	Said Reda (1969)	\$0	Ear Infection	Select	Family Dr.	Edit
	Confirmed	30-07-2022 10:40 AM to 11:04 AM	Said Reda (1969)	\$0	Eye Infection/Pinkeye	Select	Family Dr.	Edit


Showing 1 to 3 of 3 entries Previous 1 Next

This screen allows Staff/ Administrator/ Drs can view/ add/ cancel/ modify patients appointments.

Action:

- Cancel Appointment
- Un-confirm Appointment
- Upload a file related to patient
- Show payment details if any
- Process Booking if completed.

Assign or change Doctor



# Doctors Working Hours Screen

## Service Provider Working Hours

Back Save Changes

Service Provider: All Providers Branch: Riyadh

Add Services

week day < >

	Mon	Tue	Wed	Thu	Fri	Sat	Sun
8am							
8:15am							
8:30am							
8:45am							
9am	9am - 4pm Family Dr.	9am - 12:35pm Family Dr.					9am - 3pm Family Dr.
9:15am							
9:30am							
9:45am							
10am							

Select Doctor & Branch to choose, and View respective working Hours.

Drag on schedule to add, delete working schedule for selected Dr. and Branch

# Branch Blocked Time Screen

## Blocked Times

Back

Save Changes

This screen allows Administrator to create blocked times for each Branch.

10:35am - 10:45am  
cleaning works (click to remove)

11:30am - 11:45am  
luch (click to remove)

Click and Drag to create or cancel blocked time per Branch, and list the reason of this block time





# Doctors Blocked Time Screen

## Blocked Times

This screen allows Administrator to create blocked times for each Branch.

[Back](#) [Save Changes](#)

< Jul 25 – 31, 2022 >

[month](#) [week](#) [day](#)

current week

	Mon, Jul 25	Tue, Jul 26	Wed, Jul 27	Thu, Jul 28	Fri, Jul 29	Sat, Jul 30	Sun, Jul 31
all-day							
8am							
8:15am							
8:30am							
8:45am	8:45am - 9:05am preparation works (click to remove)						
9am							
9:15am							
9:30am							

Click and Drag to create or cancel blocked time Per Dr. , and list the reason of this block time



# Branches (Clinics) Screen

Search...

## Branches

This screen allows Administrator to create, edit Branches (clinics, etc...).

New Branch

Name	Description	Active	
Huron		<input checked="" type="checkbox"/>	<a href="#">Edit</a>
Baseline		<input checked="" type="checkbox"/>	<a href="#">Edit</a>
Riyadh	Riyadh Nuzha Clinic	<input checked="" type="checkbox"/>	<a href="#">Edit</a>

Shows name, and description  
Of the branch

Click to edit branch (clinic),

Click to create new branch (clinic),

Shows Status of the Branch (Active or not)

## Edit Branch Huron

[Back](#) [Save Changes](#)

Name \*

Description

Address

Phone Number

Email

Upload Picture  
 No file chosen

Interval (Minutes)

Enter Branch details

Clients Pay

Booking Payment

Active  CanBookOnline

Branch Encounter Types

AcceptVideo  AcceptPhone

ScreenCovid

Video Appt. Text

Phone Appt. Text

AcceptInPerson

In Person Appt. Text

Display Order

# Branch (Clinics ) Screen (2)


## Edit Branch Huron

Back

Save Changes


8:30am						
9am	9am - 5pm	9am - 5pm	9am - 5pm	9am - 5pm	9am - 5pm	
9:30am						
10am						10am - 2pm
10:30am						
11am						
11:30am						
12pm						
12:30pm						
1pm						
1:30pm						

Add, delete  
Branch Work Schedule



**Labs**

[+ New Lab](#)



Show 50 entries Search:

FullName	Email	Phone Number	Active	action
Test Pharamcy ( Test Pharmacy)	<a href="mailto:Pharmacy@vtest">Pharmacy@vtest</a>	789-654-1233	<input checked="" type="checkbox"/>	<a href="#">Edit</a>

Showing 1 to 1 of 1 entries

Previous 1 Next

Add, delete, activate a Pharmacy or Lab. For patient referral.

Search...


## Categories

+ New Category

### Categories



Name	Description	
General	General practionner and family Drs.	<a href="#">Edit</a>



Add, delete, edit a service category or Lab. For patient referral.

## Services

+ New Service

### Services



Category	Order	Name	Description	Service Fee (CAD)	Booking Fee (CAD)	Duration (Minutes)	Processing Break (Minutes)	Finishing Duration (Minutes)	Total Duration (Minutes)	Active	
	1	Doctor Consultation		\$0.00	\$0.00	10	0	0	10	<input checked="" type="checkbox"/>	<a href="#">Edit</a>
	1	Family Dr.	General Practitioner	\$80.00	\$20.00	10	3	2	15	<input checked="" type="checkbox"/>	<a href="#">Edit</a>
	2	ENT Dr.	ENT Dr.	\$120.00	\$30.00	15	3	2	20	<input checked="" type="checkbox"/>	<a href="#">Edit</a>
	3	Gyno. Dr.	Gynecologist	\$120.00	\$30.00	15	3	2	20	<input checked="" type="checkbox"/>	<a href="#">Edit</a>
	4	Gasto. Dr.	gastroenterologist	\$120.00	\$30.00	15	3	2	20	<input type="checkbox"/>	<a href="#">Edit</a>
	4	Gastro. Dr.	gastroenterologist	\$120.00	\$30.00	15	3	2	20	<input checked="" type="checkbox"/>	<a href="#">Edit</a>

Shows name, and description  
Cost, duration and status of each provided service

Click to edit, modify, activate or deactivate provided service

Search...

## Staff

+New Staff



Show 50 entries

Search:

FullName	Email	Phone Number	Active	action
----------	-------	--------------	--------	--------

No data available in table

Showing 0 to 0 of 0 entries

Previous Next

Shows name, contacts, status and role  
Of staff (e.g. Nurse, Administrator, etc...)



# Service Providers (Drs) Screen

## ServiceProviders

+ New Service Provider



Show 50 entries

Search:

FullName	Email	Phone Number	Active	action
Family Dr.	said1reda1@yahoo.com	519-872-8628	<input checked="" type="checkbox"/>	<a href="#">Edit</a>
General Practitioner	doctor@vtest	123-456-7893	<input checked="" type="checkbox"/>	<a href="#">Edit</a>

Showing 1 to 2 of 2 entries

Previous 1 Next

Shows name, contacts, status and role  
Of Service Providers (Drs, Specialists, Chiropractors, etc...)

Click to edit, modify, activate or deactivate  
Service Provider

# Administrator/ Dr. Clients Access Screen

Search...

## Clients



Show 50 entries

Search:

First Name	Last Name	Date of Birth	Health Card Number	Version Code	Please Select Your Status	Patient Phone Number	action
dhge	ehgehg	01-01-1990	2435-453-223	df	HC	567-654-3245	<a href="#">Edit</a>
Said	Reda	31-12-1969	1213-134-442	ee	HC	519-864-8563	<a href="#">Edit</a>

Showing 1 to 2 of 2 entries

Previous 1 Next

Shows names, contacts, status of patients

# Administrator Clients Login Screen

## Clients' Logins

+New Client



Show 50 entries

Search:

FullName	Email	Phone Number	Active	action
dhge ehgeh	<a href="mailto:ektechsolutions@gmail.com">ektechsolutions@gmail.com</a>	567-654-3245	<input type="checkbox"/>	<a href="#">Edit</a>
Said Reda	<a href="mailto:said-reda@esa-solutions.com">said-reda@esa-solutions.com</a>	519-864-8563	<input checked="" type="checkbox"/>	<a href="#">Edit</a>

Showing 1 to 2 of 2 entries

Previous 1 Next



Shows names, contacts, status of patients

[+ Upload New File](#)

**Bookings (0)**

**Files (0)**

**Payments (0)**

**Profile**

No events to display

[+ Upload New File](#)[« Back](#) 26 July 2022 [Next »](#)[Current Day](#)Show  entriesSearch: 

Name	Client	Fee	Date Uploaded	Upload By	Expiry Date	Action
No data available in table						

Showing 0 to 0 of 0 entries

[Previous](#) [Next](#)


There is no item in your cart


**Users**


Search...

 Family Dr.

 Said Reda

 General Practitioner

 dhge ehgeh

 Test Pharamcy

if you need to change your password or manage your profile, please

2022-07-26 9:02:56 PM Seen

Type a message



Payment Date From

Payment Date To

Doctor

Patient

Source

Booking Date From

Booking Date To

File Upload Date From

File Upload Date To

Booked by Pharmacy

 Search: 

client	Payment Date	Amount	Status	Ref Code	Approval Code	Card Type	Card Number	Source	Source ID
No data available in table									
<b>Total:</b>		<b>\$0</b>							

Showing 0 to 0 of 0 entries

Search... ▼

## Labs

[+New Lab](#)



Show  entries

Search:

FullName <span>↑↓</span>	Email <span>↑↓</span>	Phone Number <span>↑↓</span>	Active <span>↑↓</span>	action <span>↑↓</span>
Test Pharamcy ( Test Pharmacy)	<a href="mailto:Pharmacy@vtest">Pharmacy@vtest</a>	789-654-1233	<input checked="" type="checkbox"/>	<a href="#">Edit</a>

Showing 1 to 1 of 1 entries



## Edit Test Pharamcy

Back

Save Changes

### Account Information

First Name \*

Test

Last Name \*

Pharamcy

Email \*

Pharmacy@vtest

Phone Number \*

789-654-1233

Password \*

Repeat Password

### More Details

Pharmacy of Preference (Pharmacy Name)

Test Pharmacy

Pharmacy Address

1472 Huron Street London ON N5V 2E5

Active

[+ New Booking](#)

[« Back](#) July 26, 2022 [Next »](#)
[Current Day](#)

 Show  entries

 Search: 

↑↓	Status	Time	Name	Unpaid	Reason	Action	Assign	Edit
	<u>Validating Health Card</u>	26-07-2022 10:15 AM to 10:25 AM	<u>Said Reda (1969)</u>	\$0	Ear Infection	Select		<a href="#">Edit</a>
	<u>Confirmed</u>	26-07-2022 05:10 PM to 05:20 PM	<u>dhge ehgehg_(1990)</u>	\$0	Abdominal Pain	Select	General Practitioner	<a href="#">Edit</a>

Showing 1 to 2 of 2 entries

[Previous](#) **1** [Next](#)

Search...

New Booking

Day

« Back July 26, 2022 Next »

Current Day

Show 50 entries

Search:

↑↓	Status	Time	Name	Unpaid	Reason	Action	Assign	Edit
	<u>Validating Health Card</u>	26-07-2022 10:15 AM to 10:25 AM	<u>Said Reda (1969)</u>	\$0	Ear Infection	Select		Edit
	<u>Confirmed</u>	26-07-2022 05:10 PM to 05:20 PM	<u>dhge ehgehg (1990)</u>	\$0	Abdominal Pain	Select	General Practitioner	Edit

Showing 1 to 2 of 2 entries

Previous 1 Next

- Select
- Cancel Appt.
- Confirm
- Un-confirm Appt.
- Upload File

# Service Provider Working Hours

[Back](#)
[Save Changes](#)

Service Provider

Branch

[Add Services](#)
[week](#) [day](#) [<](#) [>](#)

	Mon	Tue	Wed	Thu	Fri	Sat	Sun
9am	9am - 5pm General Practitioner	9am - 5pm General Practitioner	9am - 5pm General Practitioner	9am - 5pm General Practitioner	9am - 5pm General Practitioner		
9:30am							
10am						10am - 2pm General Practitioner	
10:30am							
11am							
11:30am							
12pm							

# Blocked Times

[Back](#) [Save Changes](#)

Huron ▾

< Jul 25 – 31, 2022 >

[month](#) [week](#) [day](#)

current week

	Mon, Jul 25	Tue, Jul 26	Wed, Jul 27	Thu, Jul 28	Fri, Jul 29	Sat, Jul 30	Sun, Jul 31
all-day							
9am							
9:30am							
10am							
10:30am							
11am							

Search...

All Branches

All Doctors

Select Doctors

July 25 - 31, 2022

New Booking

Month

Week

3-Day


Day

Current Week

	Mon, Jul 25	Tue, Jul 26 (B:2, C:0)	Wed, Jul 27 (B:1, C:0)	Thu, Jul 28	Fri, Jul 29	Sat, Jul 30	Sun, Jul 31
8am							
8:05							
8:10							
8:15							
8:20							
8:25							
8:30							
8:35							
8:40							
8:45							
8:50							
8:55							
9am							
9:05							
9:10							
9:15							
9:20							
9:25							
9:30							
9:35							
9:40							
9:45							
9:50							
9:55							
10am							
10:05							
10:10							
10:15		10:15 Said R   V   F					
10:20							


# Patient Video Encounter Screen

**Waiting room**  
**You're the only one here**


Picture in picture 

- 
- 
- 
- 

Patient uses this screen to complete video encounter with the Doctor.



A set of navigation controls located in the bottom right corner. It features a blue directional pad with four arrows (up, down, left, right) and a central circle. Below the directional pad is a volume slider with a blue circle on the left and a plus sign on the right.

Type a message... 

- 
- 
- 

09:37



said reda


















# Patient Video Encounter Screen


**Waiting room**


You're the only one here

Picture in picture 

### Settings

Camera:

Output:  

Input:  



PTZ:  Buttons  Sliders


OK

Powered by EK Tech v1.0.0

- 
- 
- 
- 

Setting, allows to setup Audio/Video devices at patient side (Mic/Speaker/Camera).

Type a message... 



# Patient Video Encounter Screen

Waiting room

You're the only one here


**Medical Device** X

Device:

OK



Medical Video Devices (All connected video devices can be selected here, based on Dr. and Staff requests- e.g. Skin Scope, ENT Scope, etc..).




Type a message...





04:35



Doctor can view images and videos live, while patient or staff are using medical instruments (e.g. Skin camera).

- Settings icon (gear)
- Microscope icon
- U-shaped icon
- Cloud icon
- iMD icon

▶ said reda

Full-screen icon, Microphone icon, End call icon (red), Video off icon, **Share Screen** button

**Take Snapshot**

08:48



Doctor can view images and videos live, while patient or staff are using medical instruments (e.g. ENT camera).



iMD



said reda



Share Screen

Take Snapshot

# Patient Video Encounter Screen

**Waiting room**

**You're the only one here**

**Stethoscope** X

Stethoscope:

OK

- Settings (gear icon)
- Microphone (microphone icon)
- Stethoscope (stethoscope icon)
- Heart rate (heart rate icon)

Stethoscope (Connected Stethoscope can be selected here), Accordingly Dr. can listen to, and record stethoscope Audio.

Navigation controls: directional pad, central button, volume slider, and zoom in (+) button.

Control buttons: full screen, mute, and video off.

Type a message...

10:39



### Stethoscope

Stethoscope: Microphone (3- PCP-USB) ▾

Volume:

**Stop Recording**

OK

Doctor can listen to stethoscope live, while patient or staff are using medical instruments (e.g. Stethoscope).

- Settings
- Microscope
- Urgency
- Cloud
- iMD

Navigation pad with directional arrows and a central button.

Volume slider with minus and plus buttons.


said reda

Control bar with icons for full screen, microphone, call, video, and a green **Share Screen** button.

# Patient Video Encounter Screen

Waiting room

You're the only one here

Picture in picture 

## Vitals

X

Blood Pressure:  mmHg

Temperature:  C°

Pulse:  bpm

SpO2:  %

Weight:  kg

Height:  cm


Clear

Save

Vitals. Patient or on-site staff, enters vital data based on measurements then data is displayed to Doctor seamlessly.

Type a message...





Data Log. All exchanged data, including medical records are displayed in this section.

Instant chat. Patient and Doctors can have instant chats, exchange data etc...

Type a message...



Medical Camera Snapshots:

---

Stethoscope recordings:

---

Export as pdf

Export to emr

Save as zip