

MDXP-6100 Palm Patient Monitor



Patient Monitor

MDXP-6100 this monitor is capable of monitoring parameters: PR, HR(heart rate), NIBP(noninvasive blood pressure), Spo2, Resp(Respiration), TEMP(body temperature).



Color: Black and white classic color

Display: 2.8' TFT

Size: 145(H)x85(W)x37(T)

Features:

- ◆ 6 - Vital sings
 - ◆ Portable, palm size
 - ◆ One touch operation
 - ◆ APP for Android and iOS
 - ◆ Real time data transmission
 - ◆ Support RPM integration
 - ◆ Can be used while charging
 - ◆ CE0123, ISO13485 certified

Standard configuration :



1. SpO2 Sensor
2. ECG & Temp Cable
3. NIBP Cuff
4. Type-C Charging Cable
5. Charging Base