


MDXP-6750 Patient Monitor & Triage



Feature

► Six parameter module is including Spo2, Pulse Rate, ECG, NIBP and Resp functions.

It can be connected by serial RS232, Bluetooth and USB.

► Patient data upload stable, good mobility ,convenient and reliable.

Continuous monitoring facilitate customer secondary development according to specific needs ,

Accessories

1. ECG Cable
2. NIBP Cuff
3. NIBP tube
4. Temp Sensor
5. Spo2 Sensor



Features

- Spo2, Pulse Rate
- ECG, 5 leads, 7 channel
- Heart Rate
- Respiration
- Non-invasive Blood Pressure
- Temperature
- Communication: Bluetooth Dual Mode (3.0 and BLE4.0) COM, Type B

