

MDXP- GR4301 - e-stethoscope



MDXP-GR4031 is an e-stethoscope, Chestpiece with embedded amplifier, A/D converter, encoder and USB interface

The MDXP-GR4031 electronic stethoscope is designed to provide auscultatory sound via Telemedicine, computer or through mobile device platforms. This can be done in real time remotely or through sharing of recorded files.



- Uses patented piezo technology for picking up body sounds with minimal background arte- facts
- No requirements for a second doctor on transmit end. Patients can use the stethoscope on themselves
- Easy to use electronic stethoscope, no buttons or dials
- No batteries needed
- Ergonomic design