

## MDXP-9.2 Portable Ultrasound Scanner



The **USB Portable Ultrasound Scanner MDXP-9.2** provides Clear and rapid imaging for intervention guidance. Also, this USB Ultrasound Machine automatically syncs data to your private cloud for telemedicine and mobile telemedicine.

Moreover, This USB Ultrasound Scanner supports Windows Operation Systems to make its applications easy for the healthcare taker.

Last, but not least, this Portable Ultrasound Machine enables the medical to perform ultrasound scans whenever and wherever necessary...

- Small and light.
- Full digital imaging technology, crystal-clear Image.
- Windows OS Compatible.
- Abundant functions, Report Function.
- Suitable for emergencies, clinics, Outdoor Inspections And Veterinary Uses.



## Specification:

- **Scanning type:** Electronic convex array scanning
- **Scanning mode:** B, B/B, B/M, M, 4B
- **Probe type:** Convex array probe, Radius 60mm, Frequency scope: 2.0 MHz~4.5 MHz(Center Frequency 3.5MHz,
- **Elements:** 80 elements.
- **Operating System:** User's computer or tablet, Win XP, VISTA, WIN7, WIN8, Wind10.
- **Software language:** Chinese/English, switchable
- **Multi-Frequency Switch:** 2.0MHz, 2.5MHz, 3.0MHz, 3.5MHz, 4.0MHz, 4.5MHz。
- **Focus:** Electronic focus selection
- **Application Preset:** General, Gynecology, Obstetrics, Cardiology, Small parts, Urology
- **Image Processing**
- **Pre-processing:** Frame Correlation, Dynamic Rang, THI, Smoothing and sharpening
- **Post-Processing:** γ adjustment
- **Image Magnification:** 1.0, 1.2, 1.5, 2.0
- **Gain Adjustment:** Total Gain and Steps Gain adjustment
- **Display Mode:** B mode, M mode
- **Image Save Feature:** HDD Save/Preview image files (JPG)
- **Cine Loop Memory:** 512-frame memory, manual/automatic cine loop.
- **Image flip:** Left/right Mirror Horizontally Flip, Up/down Vertically Flip
- **Notation:** Date, Time, Name, Patient ID, Age, Doctor, Hospital Name; Full screen annotation edit.
- **Image report and print function.**
- **Size:** 23mm×150mm×70mm.
- **Cable Length:** 190cm.
- **Weight:** 0.2 KG.

**Certifications:** FDA CLEARED, CE, ISO13485.